



# ReadyDock:UV

## Chemical-Free disinfection

ReadyDock:UV is the worlds first product to offer chemical-free general purpose disinfection for Tablet Computers. In a medical setting, Tablets are roaming high-touch surfaces which require frequent disinfection. Chemical disinfection of a Tablet can take up to 15 minutes following a typical Hospital Disinfection SOP.



ReadyDock:UV applies germicidal light consistently to 99% of the surface area of the Motion C5 Mobile Clinical Assistant. General purpose disinfection in just 105 seconds...without chemicals.

Disinfect without the mess.



Motion Computing C5 MCA

## Features

### General Purpose Disinfection in 105 Seconds

ReadyDock:UV applies germicidal light to 99% of the surface area of your tablet for 105 seconds.

### Front Door Loading/Automated Start

Once you load your Motion C5 MCA for disinfection and close the door the disinfection cycle begins automatically. Indicators let you know when a cycle has started and is complete.

### Wall or Cart Mountable

Place ReadyDock:UV in workflow friendly locations.

### Lamp Replacement Indicators

Indicators let you know if a lamp stops functioning.

### User Replaceable Lamps

It is recommended that your lamps be replaced each year. But don't worry, you can do it yourself!

### Lamp Replacement Supplies & Service

Service agreements include automatic shipment of replacement lamps and reminders ahead of time.

**ReadyDock:UV** does not replace your standard cleaning practice, but plays a new part in methodical general purpose disinfection.

Germicidal light will not penetrate dirt and grime therefore it is highly recommended that TabletPCs be physically cleaned periodically to remove dirt and grime.

### Dimensions:

Width	14.6"	Weight	13 Lbs
Depth	3.4"	Depth with door open	11"
Height	15.9"		

### Other Specifications:

Protective circuitry prevents human exposure to germicidal light by not allowing the internal lamps to be on if the door is open.

Two 100 x 100 mm hole patterns on the rear of the device for wall or clinical cart mounting.

#### Front Panel Controls & LEDs

- Replace Left Lamp
- Replace Right Lamp
- Replace Center Internal Lamps
- UV Cycle Complete
- UV Cycle In Progress

Easy removal of germicidal lamps without disassembling ReadyDock or removing from cart or wall

Pen Holder for optional pen disinfection

Maximum Usage Cycle: 20,000 Disinfection Cycles per tablet

**For more information or to order the ReadyDock:UV contact Add-On Data**



**(978) 988-1900**  
**www.addondata.com**